

Enabling Cannulation and Smooth Stricture Passage

ABOUT ▶

DEVICE
FEATURES ▶

IMPORTANT
PROCEDURAL
NOTE ▶

THERAPEUTIC
ACCESSORIES ▶

ORDERING ▶

*Optimizing wireguided cannulation for the
treatment of pancreatico-biliary diseases*

NaviPro™ Hydrophilic Guidewire

ABOUT ▶

DEVICE
FEATURES ▶

IMPORTANT
PROCEDURAL
NOTE ▶

THERAPEUTIC
ACCESSORIES ▶

ORDERING ▶



About the NaviPro Guidewire

The NaviPro Guidewire is a 100 percent Hydrophilic guidewire designed for cannulation and smooth passage through tortuous anatomy.

This family of performance guidewires is designed for physicians who prefer using over-the-wire access and therapeutic devices for endoscopic retrograde cholangio-pancreatography (ERCP).

“Guidewires have become an integral part of ERCP. A hydrophilic wire, particularly the tip, is ideal for gaining access into different parts of the ductal systems. For ninety percent of our cases we use a hydrophilic wire for access, cannulation and placement of therapeutic devices.”

— **Dr. David Carr-Locke**
Beth Israel Medical Center, New York, NY

Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT

DEVICE
FEATURES

IMPORTANT
PROCEDURAL
NOTE

THERAPEUTIC
ACCESSORIES

ORDERING



1 of 2

Device Features

NaviPro Hydrophilic Guidewires are available in two tip configurations to facilitate cannulation of the biliary ducts in a variety of therapeutic cases.

ANGLED

STRAIGHT

Each guidewire has a working length of 260cm and is available in three outer diameter sizes:

- .018" / 0.46mm
- .025" / 0.64mm
- .035" / 0.89mm

Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT

DEVICE
FEATURES

IMPORTANT
PROCEDURAL
NOTE

THERAPEUTIC
ACCESSORIES

ORDERING



◀ 2 of 2

Device Features

NaviPro Hydrophilic Guidewire offers a one-to-one torque response for selective cannulation. The highly radiopaque guidewire provides excellent visibility to facilitate ductal location and accurate positioning during ERCP.

The torqueing device aids the physician or technician providing a one-to-one torque response.

**Torqueing device sold separately.*



Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT

DEVICE
FEATURES

IMPORTANT
PROCEDURAL
NOTE

THERAPEUTIC
ACCESSORIES

ORDERING



Important Procedural Note

To activate the Hydrophilic coating, prior to removing the guidewire from its dispenser, inject a sterile saline solution into the luer lock hub using a syringe.



Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT



DEVICE
FEATURES



IMPORTANT
PROCEDURAL
NOTE



THERAPEUTIC
ACCESSORIES



ORDERING



1 of 2

Therapeutic Accessories

Boston Scientific's wide range of endoscopic access and therapeutic devices provides physicians with options to implement the most effective course of treatment.

Through innovation and continuous improvement, Boston Scientific offers a range of therapeutic accessories for the diagnosis and treatment of pancreatico-biliary stone and stricture management.



Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT

DEVICE
FEATURES

IMPORTANT
PROCEDURAL
NOTE

THERAPEUTIC
ACCESSORIES


ORDERING




◀ 2 of 2

Therapeutic Accessories


The NaviPro Hydrophilic Guidewire works effectively with these devices:




Flexima™ Biliary Stent




Contour™ Cannula
Fluoro Tip™ Cannula
Tandem™ XL Triple-Lumen Cannula



CRE™ Balloon for Dilation
Assisted Stone Extraction



Extractor™ Pro Retrieval
Balloon XL and DL



TRUEtome™ Cannulating
Sphincterotome

Visit: www.bostonscientific.com/endo-resources for product information and links to procedural videos.

Boston
Scientific

NaviPro™ Hydrophilic Guidewire

ABOUT

DEVICE
FEATURES

IMPORTANT
PROCEDURAL
NOTE

THERAPEUTIC
ACCESSORIES

ORDERING



Ordering Information

The NaviPro Hydrophilic Guidewire length is 260cm.

Product Code	Guidewire Diameter (inches / mm)	Tip	Packaging
M005 5619 1	.018 / 0.46	Angled	Box 5
M005 5618 1	.018 / 0.46	Straight	Box 5
M005 5621 1	.025 / 0.64	Angled	Box 5
M005 5620 1	.025 / 0.64	Straight	Box 5
M005 5623 1	.035 / 0.89	Angled	Box 5
M005 5622 1	.035 / 0.89	Straight	Box 5
M005 5625 1	.035 / 0.89	Angled Stiff	Box 5
M005 5624 1	.035 / 0.89	Straight Stiff	Box 5

**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
One Boston Scientific Place
Natick, MA 01760-1537
www.bostonscientific.com/endo-resources

Ordering Information
1.888.272.1001

©2013 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-151302-AB May 2013