



One Less Worry.



Orca™ Sterile Single Use
Air/Water and Suction Valves

Seal™ Biopsy Valve &
Hydra™ Disposable Auxiliary
Water Jet Connector

Ordering
Information

Orca™

Sterile Single Use Air/Water and Suction Valves



Enclosed spring

Designed to prevent entrapment of fluid or debris which could compromise function and cause valve "sticking"

Silicone free

In accordance with manufacturer's recommendations

Recessed, matte finish finger grip

Helps provide greater gripping capability and minimizes finger slippage

Precision manufactured shaft and symmetrical center hole

Consistent fit, function and performance



Snap-fit boot

Designed to secure attachment of valve to cylinder to help prevent leakage



Large, Flat, Color-Coded Finger Perch

Provides comfort and easy valve identification

Orca™
Sterile Single Use
Air/Water and Suction Valves

Seal™ Biopsy Valve &
Hydra™ Water Jet Connector

Competitive
Comparisons

Ordering
Information

Seal™


Single Use Biopsy Valve

- Single Use to reduce the risk of cross-contamination between patients
- Designed to minimize leakage and maintains scope pressure while passing device

Hydra™

Disposable Auxiliary Water Jet Connector

- Single Use - no lost components to replace
- One-way Valve designed to reduce the risk of cross-contamination
- Metal Threading designed to minimize potential to cross thread when connecting to auxiliary water jet




Orca™
Sterile Single Use
Air/Water and Suction Valves

Seal™ Biopsy Valve &
Hydra™ Water Jet Connector

Competitive
Comparisons

Ordering
Information

Ordering Information



Order Number	Product Description	Orca™ Valves	Seal™ Biopsy Valve	Hydra™ Disposable Auxiliary Water Jet Connector	Scope Compatibility	Units
SUV-617-50	Orca Sterile, Single Use Air/Water and Suction Valves	✓			Olympus™ 160/180/190 Series	Box 50
SUV-628-50	OrcaPod: Orca Sterile, Single Use Air/Water and Suction Valves with Seal™ Biopsy Valve	✓	✓		Olympus 160/180/190 Series	Box 50
SUV-629-50	OrcaPod: Orca Sterile, Single Use Air/Water and Suction Valves, Seal Biopsy Valve and Hydra Water Jet Connector	✓	✓	✓	Olympus 160/180/190 Series with auxiliary water port	Box 50
SBC-501-100	Seal Biopsy Valve for Olympus & Fuji™		✓		Olympus 160/180/190 and Fuji 530/570/600 and G5	Box 100
SIT-361-100	Hydra Disposable Auxiliary Water Jet Connector			✓	Olympus 160/180/190	Box 100

Sterile Single Use Air/Water and Suction Valves

Seal™ Biopsy Valve & Hydra™ Water Jet Connector

Competitive Comparisons

Ordering Information



**Boston
Scientific**
Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

©2018 Boston Scientific Corporation
or its affiliates. All rights reserved.

ENDO-470733-AB

All photographs are owned by Boston Scientific. All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

Rx Only.

CAUTION: The law restricts this device to sale by or on the order of a physician. Information for use only in countries with applicable health authority registrations.

Materials not intended for use in France.