

Polypectomy Snares

**Boston
Scientific**
Advancing science for life™



Polypectomy Snares

Available in a broad range of loop sizes and shape configurations.

Ordering Information

CAPTIVATOR™ II SINGLE-USE SNARE

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561291	Extra Large Rounded	Stiff	33	2.4	240	2.8	Box 10

CAPTIVATOR™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562301	Small Oval	Stiff	13	2.4	240	2.8	Box 10
M00562321	Medium Oval	Stiff	27	2.4	240	2.8	Box 10
M00561311	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 10
M00562391	Large Oval (thin wire)	Flexible	30	2.4	240	2.8	Box 10
M00562451	Small Hexagonal	Stiff	13	2.4	240	2.8	Box 10
M00562341	Medium Hexagonal	Stiff	27	2.4	240	2.8	Box 10
M00562371	Medium Crescent	Stiff	27	2.4	240	2.8	Box 10

SINGLE-USE ROTATABLE SNARES*

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00561821	Small Oval	Medium Stiff	13	2.4	240	2.8	Box 10
M00561831	Medium Oval	Medium Stiff	20	2.4	240	2.8	Box 10

* With Star Catheter Technology

SENSATION™ SINGLE-USE SHORT THROW SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562691	Small Oval	Flexible	13	2.4	240	2.8	Box 5
M00562692	Small Oval	Flexible	13	2.4	240	2.8	Box 20
M00562693	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00560321	Small Oval	Medium Stiff	13	2.4	240	2.8	Box 5
M00562671	Medium Oval	Flexible	27	2.4	240	2.8	Box 5
M00562672	Medium Oval	Flexible	27	2.4	240	2.8	Box 20
M00562673	Medium Oval	Flexible	27	2.4	240	2.8	Box 40
M00560311	Medium Oval	Medium Stiff	27	2.4	240	2.8	Box 5
M00562651	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 5
M00562652	Large Oval	Medium Stiff	30	2.4	240	2.8	Box 20
M00560331	Medium Crescent	Medium Stiff	27	2.4	240	2.8	Box 5

CAPTIFLEX™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562471	Extra Small Oval	Flexible	11	2.4	240	2.8	Box 10
M00562421	Small Oval	Flexible	13	2.4	240	2.8	Box 10
M00562422	Small Oval	Flexible	13	2.4	240	2.8	Box 40
M00562401	Medium Oval	Flexible	27	2.4	240	2.8	Box 10
M00562402	Medium Oval	Flexible	27	2.4	240	2.8	Box 40

PROFILE™ SINGLE-USE SNARES

Order Number	Shape	Stiffness	Loop Width (mm)	Sheath O.D. (mm)	Working Length (cm)	Min. Working Channel (mm)	Units
M00562531	Extra Small Oval	Flexible	11	1.9	240	2.0	Box 10
M00562551	Small Oval	Flexible	13	1.9	240	2.0	Box 10
M00562571	Medium Oval	Flexible	27	1.9	240	2.0	Box 10

Boston Scientific

Advancing science for life™

Boston Scientific Corporation
300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com/gastro

Ordering Information
1.888.272.1001

©2016 Boston Scientific Corporation or its affiliates. All rights reserved.

ENDO-145002-AB August 2016

All photographs taken by Boston Scientific.

All trademarks are the property of their respective owners.

Indications, Contraindications, Warnings and Instructions for Use can be found in the product labeling supplied with each device.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION: The law restricts these devices to sale by or on the order of a physician. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.