


PERIVAC™

PERICARDIOCENTESIS KITS

PRE-ASSEMBLED KIT FOR PERICARDIAL DRAINAGE PROCEDURES



PROVIDES THE ESSENTIALS FOR PERICARDIAL DRAINAGE/ASPIRATION

PERIVAC™

PERICARDIOCENTESIS KITS

Ordering Information

Model	Description	Quantity
M004 PERIVAC1P 0	PERIVAC Kit with 8.3F Pigtail Catheter	1 Kit per box
M004 PERIVAC1S 0	PERIVAC Kit with 8.3F Straight Catheter	1 Kit per box
M004 400 1	Spare 1000cc Drainage Bag (individually sterilized)	3 Units per box

To Order, Call Boston Scientific
Customer Service at 1-888-272-1001

PERIVAC™ KIT from Boston Scientific

INDICATIONS FOR USE: The Pericardiocentesis Kit is intended to withdraw fluid from the pericardial sac (membrane surrounding the heart) to evaluate the fluid for diagnostic purposes or to remove excess fluid causing compression of the heart.

CONTRAINDICATIONS: There are no known contraindications for pericardial aspiration, although recurrent effusion or unresolved cardiac tamponade may warrant surgical intervention.

COMPLICATIONS/RISKS: Complications and/or risks of pericardiocentesis procedures can include catheter occlusion, infection, cardiac tamponade from puncture, atrial and ventricular arrhythmias, laceration of coronary vessel or chamber wall, air injection into cardiac chambers, pneumothorax, hemothorax, and hypotension (perhaps reflexogenic). Perforation of the stomach, colon, right atrium and/or ventricle is also possible.

WARNINGS: Rx ONLY Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

- This device is intended for single use only.
- Store in a cool dry place.
- A guide wire should never be advanced or withdrawn against resistance. The cause of the resistance should be determined before continuing.

PRECAUTIONS: Pericardiocentesis should only be performed by clinicians well-versed in the technique and aware of its risks and benefits.

EP-992505-AA

Boston Scientific is a distributor of the PERIVAC kit, manufactured by Merit Medical.

All trademarks are the property of their respective owners.

**Boston
Scientific**

Advancing science for life™

Rhythm Management

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

Medical Professionals:
1.800.CARDIAC (227.3422)
Customer Service:
1-888-272-1001

© 2021 Boston Scientific Corporation
or its affiliates. All rights reserved.

EP-143805-AC