

ACCUSTICK™ II Introducer Systems With Radiopaque Marker

FOR THE INTRODUCTION OF .038" GUIDEWIRES IN NON-VASCULAR PROCEDURES

RADIOPAQUE MARKER

RO marker positioned at tip of sheath to enhance visualization for more precise placement

LARGE SHEATH END HOLE INNER DIAMETER

Facilitates placement of .038" working guidewire alongside .018" guidewire

CONVENIENT, INTEGRATED SYSTEM

- 21 Gauge diagnostic needle designed to reduce tissue trauma
- Option of .018" kink-resistant nitinol guidewire is designed to provide strength during tough procedures
- Coaxial sheath/dilator assembly with locking stiffening cannula is designed for over-the-wire placement
- .038" heavy duty guidewire in a choice of "J" tip or straight tip provides options as a working wire

SHEATH/DILATOR ASSEMBLY PROXIMAL TAPER

For ease of introduction and placement

DILATOR REFERENCE MARK

Facilitates alignment of sheath and dilator tips for accurate wire placement



ACCUSTICK™ II Introducer Systems With Radiopaque Marker

UPN	ORDER NUMBER	DESCRIPTION	QUANTITY
M00127020*	20-702	AccuStick II Introducer System with -1 Extra Stiff, Stainless Steel, 60 cm, .018" Guidewire with 7.5 cm Floppy Platinum Tip -1 Heavy Duty, Stainless Steel, 145 cm, .038" Guidewire with 3mm "J" Tip	1
M001207030*	20-703	AccuStick II Introducer System with -1 Extra Stiff, Stainless Steel, 60 cm, .018" Guidewire with 7.5 cm Floppy Platinum Tip -1 Heavy Duty, Stainless Steel, 145 cm, .038" Guidewire with Straight Tip	1
M001207100*	20-710	AccuStick II Introducer System with -1 Nitinol, 60 cm, .018" Guidewire with 6 cm Floppy Platinum Tip -1 Heavy Duty, Stainless Steel, 145 cm, .038" Guidewire with 3 mm "J" Tip	1
M001207040*	20-704	AccuStick II Introducer System with -1 Extra Stiff, Stainless Steel, 60 cm, .018" Guidewire with 7.5 cm Floppy Platinum Tip	1
M001207050*	20-705	AccuStick II Introducer System without Guidewires	1
M001201081	20-108	Stainless Steel, 60 cm, .018" Guidewire with 7.5 cm Floppy Platinum Tip	5
M001207111	20-711	Nitinol, 60 cm, .018" Guidewire with 6 cm Floppy Platinum Tip	5
M001206041	20-604	21 Gauge 15 cm Diagnostic Needle with Stylet	10

*Each system includes: One AccuStick II Introducer (Coaxial 6F Sheath with RO Marker/4F Dilator. 20 cm) with Locking Stiffening Cannula; One 21 Gauge, 15 cm Diagnostic Needle with Stylet.



ACCUSTICK™ II INTRODUCER SYSTEM

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INTENDED USE/INDICATIONS FOR USE: The AccuStick™ II Introducer System with radiopaque marker facilitates introduction and placement of a guidewire.

CONTRAINDICATIONS: None known. **PRECAUTIONS:** None known. **WARNINGS:** None known. **ADVERSE EVENTS:** Potential risks exist for serious complications to include: • Perforation of a vessel or viscus • Laceration of a vessel or viscus • Bleeding • Wire or catheter embolism • Extravasation • Hematoma • Hemothorax • Hydrothorax • Inflammation, necrosis or scarring • Risks normally associated with percutaneous interventional procedures • Pain in region • Skin infection • Edema
These and other complications are well documented in medical literature. Use of the AccuStick II Introducer System with radiopaque. **90960661 Rev/Ver. AB**

Peripheral Interventions

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information
contact customer service at 1.888.272.1001.

© 2020 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-390501-AB